















Glossary of Symbols

Symbol	Title of Symbol	Description of Symbol	Standard or Reference
	Manufacturer	Indicates the medical device manufacturer	ISO 15223-1:2021 Clause 5.1.1
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified	ISO 15223-1:2021 Clause 5.1.6
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified	ISO 15223-1:2021 Clause 5.1.5
	Material	Indicates material composition of implant	n/a
	Use-by date (YYYY-MM-DD)	Indicates the date after which the medical device is not to be used	ISO 15223-1:2021 Clause 5.1.4
	Date of manufacture (YYYY-MM-DD)	Indicates the date when the medical device was manufactured (package date)	ISO 15223-1:2021 Clause 5.1.3
	Prescription Use Only	Requires prescription in the United States	21 CFR 801.109
	Sterilized using irradiation	Indicates a medical device that has been sterilized using irradiation	ISO 15223-1:2021 Clause 5.2.4
	Do not re-sterilize	Indicates a medical device that is not to be re-sterilized	ISO 15223-1:2021 Clause 5.2.6
	Single Sterile Barrier System with Protective Packaging Inside	Indicates a single sterile barrier system with protective packaging inside.	ISO 15223-1:2021 Clause 5.2.13
	Do not use if package is damaged and consult instructions for use	Indicates that a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information	ISO 15223-1:2021 Clause 5.2.8
	Do not re-use	Indicates a medical device that is intended for one single use only	ISO 15223-1:2021 Clause 5.4.2
	Consult instructions for use	Indicates the need for the user to consult the instructions for use	ISO 15223-1:2021 Clause 5.4.3
	Caution: Read all warnings and precautions in instructions for use	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself	ISO 15223-1:2021 Clause 5.4.4